


Anti-Gabarap Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O95166

Item Number: CM0016952



Introduction

Rabbit polyclonal antibody to GABARAP, validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, and rat. Ideal for autophagy and GABA receptor trafficking research.

[Learn More](#)

Parameter	Value
Product Name	GABARAP Rabbit pAb
Remarks/Alias	MM46; ATG8A; GABARAP-a; GABARAP
Species	Human
Gene ID (Human)	11337
Gene ID	11337
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human GABARAP (NP_009209.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95166
Protein Weight	14kDa
Shipping	Ice bag