


Anti-Gabarapl1 Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9H0R8

Item Number: CM0005460



Introduction

Recombinant rabbit monoclonal antibody to human GABARAPL1, validated for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat. Targets a protein involved in autophagy and receptor trafficking. Ideal for neuroscience and cell biology studies.

[Learn More](#)

Parameter	Value
Product Name	GABARAPL1 Rabbit mAb
Remarks/Alias	GABARAPL1; APG8-LIKE; APG8L; ATG8; ATG8B; ATG8L; GEC1; GABA type A receptor associated protein like 1
Species	Human
Gene ID	23710
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human GABARAPL1. (NP_113600.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H0R8
Protein Weight	14kDa/17kDa
Shipping	Ice bag