


## Anti-Gabarapl1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9H0R8

Item Number: CM0028183



### Introduction

Rabbit polyclonal antibody targeting GABARAPL1 (GEC1, ATG8L). Validated for Western blot and ELISA; cross-reacts with human, mouse, and rat. Investigates autophagy and opioid receptor trafficking.

[Learn More](#)

| Parameter        | Specification                                                                                                                  |
|------------------|--------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | GABARAPL1 Rabbit pAb                                                                                                           |
| Remarks/Alias    | ATG8, GEC1, APG8L, ATG8B, ATG8L, APG8-LIKE, GABARAPL1                                                                          |
| Species          | Human                                                                                                                          |
| GeneID (Human)   | 23710                                                                                                                          |
| GeneID           | 23710                                                                                                                          |
| Immunogen        | Synthetic peptide<br>A synthetic peptide corresponding to a sequence within amino acids 1-50 of human GABARAPL1 (NP_113600.1). |
| Source           | Rabbit                                                                                                                         |
| Category         | Polyclonal Antibodies                                                                                                          |
| Application      | WB, ELISA                                                                                                                      |
| Cross Reactivity | Human, Mouse, Rat                                                                                                              |
| SWISS            | Q9H0R8                                                                                                                         |
| Protein Weight   | 14kDa                                                                                                                          |
| Shipping         | Ice bag                                                                                                                        |