


Anti-Gabarapl2 Polyclonal Antibody For Wb, Elisa - P60520

Item Number: CM0020665



Introduction

Rabbit polyclonal antibody against human GABARAPL2, suitable for WB and ELISA. Reacts with human and mouse samples. Targets a protein involved in autophagy and intra-Golgi traffic; ideal for cell biology research.

[Learn More](#)

Parameter	Value
Product Name	GABARAPL2 Rabbit pAb
Remarks/Alias	ATG8; GEF2; ATG8C; GEF-2; GATE16; GATE-16; GABARAPL2
Species	Human
Gene ID	11345
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-50 of human GABARAPL2 (NP_009216.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P60520
Protein Weight	14kDa
Shipping	Ice bag