

Anti-Gabarapl2 Rabbit Monoclonal Antibody For Wb, Elisa - P60520

Item Number: CM0010528



Introduction

Rabbit monoclonal antibody targeting human GABARAPL2 (GATE-16). Validated for WB and ELISA; cross-reacts with human, mouse, rat. 14 kDa, immunogen aa 38-117. Suitable for autophagy and intra-Golgi transport studies.

[Learn More](#)

Parameter	Value
Product Name	GABARAPL2 Rabbit mAb
Remarks/Alias	ATG8, GEF2, ATG8C, GEF-2, GATE16, GATE-16, FLC3A, GABARAPL2
Species	Human
GeneID (Human)	11345
GeneID	11345
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 38-117 of human GABARAPL2 (P60520).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P60520
Protein Weight	14kDa
Shipping	Ice bag