

Anti-Gabpa Polyclonal Antibody For Wb, Elisa - Q06546

Item Number: CM0017056

Introduction

Rabbit polyclonal antibody targeting human GA-binding protein alpha chain (GABPA/NRF2 α). Suitable for Western blot and ELISA, with cross-reactivity to human and rat. Ideal for transcriptional regulation research.

[Learn More](#)

Parameter	Value
Product Name	GABPA Rabbit pAb
Remarks/Alias	NFT2; NRF2; NRF2A; E4TF1A; E4TF1-60; RCH04A07; GABPA
Species	Human
Gene ID (Human)	2551
Gene ID	2551
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 50-150 of human GABPA (NP_002031.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q06546
Protein Weight	51kDa
Shipping	Ice bag