

Anti-Gabpa Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q06546

Item Number: CM0027767

Introduction

Rabbit polyclonal antibody against human GABPA (GA-binding protein alpha chain), validated for WB, IHC-P, IF/ICC, and ELISA applications. Cross-reacts with human, mouse, and rat samples. Research use only.

[Learn More](#)

Parameter	Value
Product Name	GABPA Rabbit pAb
Remarks/Alias	NFT2; NRF2; NRF2A; E4TF1A; E4TF1-60; RCH04A07; GABPA
Species	Human
Gene ID (Human)	2551
Gene ID	2551
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human GABPA (NP_002031.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q06546
Protein Weight	51kDa
Shipping	Ice bag