

Anti-Gabpb1 Polyclonal Antibody For Wb And Elisa - Q06547

Item Number: CM0029636

Introduction

Rabbit polyclonal antibody to GABPB1, a transcription factor regulating nuclear-encoded mitochondrial genes. Validated for WB and ELISA; cross-reacts with human and mouse; 42kDa.

[Learn More](#)

Parameter	Value
Product Name	GABPB1 Rabbit pAb
Remarks/Alias	E4TF1; GABPB; BABPB2; E4TF1B; GABPB2; NRF2B1; NRF2B2; GABPB-1; E4TF1-47; E4TF1-53; GABPB1
Species	Human
GeneID (Human)	2553
GeneID	2553
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 231-395 of human GABPB1 (NP_005245.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q06547
Protein Weight	42kDa
Shipping	Ice bag