

Anti-Gabpb2 Polyclonal Antibody For Wb, Elisa - Q8Tak5

Item Number: CM0030930

Introduction

Rabbit polyclonal antibody against human GABPB2, validated for WB and ELISA. Applicable to mouse and rat samples. Ideal for studying transcription regulation involving GA-rich repeats.

[Learn More](#)

Parameter	Value
Product Name	GABPB2 Rabbit pAb
Remarks/Alias	GABPB-2; GABPB2
Species	Human
Gene ID (Human)	126626
Gene ID	126626
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 239-448 of human GABPB2 (NP_653219.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q8TAK5
Protein Weight	49kDa
Shipping	Ice bag