

Anti-Gabra1 Polyclonal Antibody For Wb, Elisa - P14867

Item Number: CM0014013



Introduction

High-quality rabbit polyclonal antibody against human GABRA1 (GABAAR alpha-1). Suitable for WB and ELISA with cross-reactivity to mouse and rat. Ideal for inhibitory neurotransmission and epilepsy research.

[Learn More](#)

Parameter	Value
Product Name	GABRA1 Rabbit pAb
Remarks/Alias	EJM; ECA4; EJM5; DEE19; EIEE19; GABRA1
Species	Human
Gene ID (Human)	2554
Gene ID	2554
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 28-251 of human GABRA1 (NP_000797.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P14867
Protein Weight	52 kDa
Shipping	Ice bag