

Anti-Gabra2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P47869

Item Number: CM0022610



Introduction

Rabbit polyclonal antibody against the GABA(A) receptor alpha-2 subunit, validated for Western blot, IHC-P, and ELISA in human and mouse samples. Ideal for neuroscience research and epilepsy studies.

[Learn More](#)

Parameter	Value
Product Name	GABRA2 Rabbit pAb
Remarks/Alias	DEE78; EIEE78; GABRA2
Species	Human
Gene ID	2555 (Human)
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 300-400 of human GABRA2 (NP_000798.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P47869
Protein Weight	51kDa
Shipping	Ice bag