


Anti-Gabra3 Rabbit Polyclonal Antibody For Wb And Elisa - P34903

Item Number: CM0015917



Introduction

Anti-GABRA3 rabbit polyclonal antibody validated for WB and ELISA. Reacts with human, mouse, and rat GABRA3. Useful for inhibitory neurotransmitter and ion channel research.

[Learn More](#)

Parameter	Value
Product Name	GABRA3 Rabbit pAb
Remarks/Alias	EPILX2; GABRA3
Species	Human
GeneID (Human)	2556
GeneID	2556
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-276 of human GABRA3 (NP_000799.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P34903
Protein Weight	55kDa
Shipping	Ice bag