

Anti-Gabra5 Polyclonal Antibody For Wb, Elisa - P31644

Item Number: CM0017253



Introduction

Rabbit polyclonal antibody targeting the gamma-aminobutyric acid receptor subunit alpha-5 (GABRA5). Suitable for WB and ELISA with reactivity across human, mouse, and rat. A valuable tool for investigating GABAergic signaling and neurological disorders.

[Learn More](#)

Parameter	Value
Product Name	GABRA5 Rabbit pAb
Remarks/Alias	DEE79; EIEE79; GABRA5
Species	Human
GeneID (Human)	2558
GeneID	2558
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 32-259 of human GABRA5 (NP_000801.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31644
Protein Weight	52kDa
Shipping	Ice bag