

Anti-Gabra5 Rabbit Polyclonal Antibody For Wb And Elisa - P31644

Item Number: CM0030552



Introduction

High-quality rabbit polyclonal antibody against GABA(A) receptor alpha-5 subunit (GABRA5). Validated for Western blot and ELISA applications, with cross-reactivity to mouse and rat. Ideal for neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	GABRA5 Rabbit pAb
Remarks/Alias	DEE79; EIEE79; GABRA5
Species	Human
Gene ID (Human)	2558
Gene ID	2558
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P31644
Protein Weight	52kDa
Shipping	Ice bag