

Anti-Gabrd Rabbit Polyclonal Antibody For Wb, Elisa - O14764

Item Number: CM0024062



Introduction

Rabbit polyclonal antibody targeting human GABRD (GABA(A) receptor delta subunit), validated for WB and ELISA. Cross-reacts with rat. Ideal for research on tonic GABAergic inhibition, sleep regulation, and epilepsy.

[Learn More](#)

Parameter	Value
Product Name	GABRD Rabbit pAb
Remarks/Alias	EJM7; EIG10; GEFSP5; GABRD
Species	Human
GeneID (human)	2563
GeneID	2563
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 31-248 of human GABRD (NP_000806.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	O14764
Protein Weight	51kDa
Shipping	Ice bag