

Anti-Gabrg1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8N1C3

Item Number: CM0015193

Introduction

Rabbit polyclonal antibody targeting human GABRG1 (GABA(A) receptor gamma-1 subunit). Suitable for WB, IF/ICC, ELISA; cross-reacts with human, mouse, rat. Ideal for inhibitory neurotransmitter receptor studies.

[Learn More](#)

Parameter	Value
Product Name	GABRG1 Rabbit pAb
Remarks/Alias	GABRG1
Species	Human
GeneID (Human)	2565
GeneID	2565
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 200-260 of human GABRG1 (NP_775807.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N1C3
Protein Weight	54kDa
Shipping	Ice bag