


Anti-Gabrr1 Rabbit Polyclonal Antibody For Wb And Elisa - P24046

Item Number: CM0026692



Introduction

Anti-GABRR1 rabbit polyclonal antibody validated for WB and ELISA. Recognizes human, mouse, rat GABRR1 (rho-1 GABA receptor subunit). Suitable for retinal and cerebellar neurotransmission research.

[Learn More](#)

Parameter	Value
Product Name	GABRR1 Rabbit pAb
Remarks/Alias	GABRR1
Species	Human
GeneID (Human)	2569
GeneID	2569
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 365-465 of human GABRR1 (NP_002033.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P24046
Protein Weight	56kDa
Shipping	Ice bag