


Anti-Gad65/Gad2 Monoclonal Antibody For Wb, Ip, Elisa - Q05329

Item Number: CM0006565



Introduction

Recombinant rabbit monoclonal antibody targeting human GAD65/GAD2 (glutamate decarboxylase 2). Validated for WB, IP, ELISA with mouse and rat cross-reactivity. Ideal for GABA biosynthesis and neurological research.

[Learn More](#)

Parameter	Value
Product Name	GAD65/GAD2 Rabbit mAb
Remarks/Alias	GAD65; GAD65/GAD2
Species	Human
Gene ID (Human)	2572
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-125 of human GAD65/GAD2 (NP_000809.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q05329
Protein Weight	65kDa
Shipping	Ice bag