


Anti-Gad65/Gad2 Mouse Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q05329

Item Number: CM0009906



Introduction

Mouse monoclonal antibody against human GAD65/GAD2, validated for WB, IF-P, IHC-P, ELISA. Cross-reacts with mouse and rat. Targets the 65kDa glutamate decarboxylase isoform involved in GABA biosynthesis.

[Learn More](#)

Parameter	Value
Product Name	GAD65/GAD2 Mouse mAb
Remarks/Alias	GAD65
Species	Human
GeneID (Human)	2572
GeneID	2572
Immunogen	Recombinant protein of human GAD65/GAD2
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q05329
Protein Weight	65kDa
Shipping	Ice bag