


Anti-Gad65/Gad2 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - Q05329

Item Number: CM0015687



Introduction

High-quality rabbit polyclonal antibody against human GAD65/GAD2, a GABA-synthesizing enzyme. Validated for WB, IF-P, and ELISA with cross-reactivity to mouse and rat. Suitable for neuroscience and autoimmunity studies.

[Learn More](#)

Parameter	Value
Product Name	GAD65/GAD2 Rabbit pAb
Remarks/Alias	GAD65; GAD65/GAD2
Species	Human
Gene ID (Human)	2572
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 401-585 of human GAD65/GAD2 (NP_001127838.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q05329
Protein Weight	65 kDa
Shipping	Ice bag