


Anti-Gad67/Gad1 Polyclonal Antibody For Wb, If-P, Elisa - Q99259

Item Number: CM0030383



Introduction

High-quality rabbit polyclonal antibody targeting human GAD67/GAD1 (GABA synthesis enzyme) for WB, IF-P, ELISA. Reacts with human, mouse, rat. Ideal for neuroscience and diabetes research.

[Learn More](#)

Parameter	Value
Product Name	GAD67/GAD1 Rabbit pAb
Remarks/Alias	GAD; SCP; CPSQ1; DEE89
Species	Human
GeneID (Human)	2571
GeneID	2571
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 141-240 of human GAD67/GAD1 (NP_000808.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99259
Protein Weight	67kDa
Shipping	Ice bag