


Anti-Gad67/Gad1 Rabbit Monoclonal Antibody For Wb And Elisa - Q99259

Item Number: CM0003224



Introduction

High-quality anti-GAD67/GAD1 rabbit monoclonal antibody validated for Western blot and ELISA. Detects human, mouse, and rat GAD67. Ideal for neuroscience, neurobiology, and GABA metabolism research.

[Learn More](#)

Parameter	Value
Product Name	GAD67/GAD1 Rabbit mAb
Remarks/Alias	GAD; SCP; CPSQ1; DEE89; GAD67/GAD1
Species	Human
GeneID (Human)	2571
GeneID	2571
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human GAD67/GAD1 (Q99259).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99259
Protein Weight	67kDa
Shipping	Ice bag