


Anti-Gadd45G Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O95257

Item Number: CM0014724



Introduction

Rabbit polyclonal antibody against GADD45G (CR6/DDIT2), validated for WB, IHC-P, ELISA in human, mouse, rat. Involved in growth arrest, apoptosis, and stress-responsive MAPK signaling.

[Learn More](#)

Parameter	Specification
Product Name	GADD45G Rabbit pAb
Remarks/Alias	CR6; DDIT2; GRP17; GADD45gamma; GADD45G
Species	Human
Gene ID (Human)	10912
Gene ID	10912
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 90-159 of human GADD45G (NP_006696.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95257
Protein Weight	17kDa
Shipping	Ice bag