

Anti-Galactosidase Alpha (Gla) Monoclonal Antibody For Wb, Elisa - P06280

Item Number: CM0008550



Introduction

Recombinant rabbit monoclonal antibody against human Galactosidase alpha (GLA), suitable for WB and ELISA. Cross-reacts with mouse and rat. Ideal for glycosphingolipid metabolism and Fabry disease research.

[Learn More](#)

Parameter	Value
Product Name	Galactosidase alpha (GLA) Rabbit mAb
Remarks/Alias	GALA; Galactosidase alpha (GLA)
Species	Human
GeneID (Human)	2717
GeneID	2717
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 50-150 of human Galactosidase alpha (GLA) (GLA) (P06280).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06280
Protein Weight	49kDa
Shipping	Ice bag