

Anti-Galactosidase Alpha (Gla) Polyclonal Antibody For Wb, Elisa - P06280

Item Number: CM0018286



Introduction

Rabbit polyclonal antibody targeting human alpha-galactosidase A (GLA), suitable for Western blot and ELISA. Supports research on lysosomal glycosphingolipid metabolism and Fabry disease.

[Learn More](#)

Parameter	Value
Product Name	Galactosidase alpha (GLA) Rabbit pAb
Remarks/Alias	GALA; Galactosidase alpha (GLA)
Species	Human
GeneID (Human)	2717
GeneID	2717
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 150-429 of human Galactosidase alpha (GLA) (NP_000160.1))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P06280
Protein Weight	49kDa
Shipping	Ice bag