

Anti-Galactosidase Alpha (Gla) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P06280

Item Number: CM0023347



Introduction

Rabbit polyclonal antibody against human Galactosidase alpha (GLA), validated for WB, IF/ICC, and ELISA. Cross-reacts with human and mouse. Protein weight 49kDa.

[Learn More](#)

Parameter	Value
Product Name	Galactosidase alpha (GLA) Rabbit pAb
Remarks/Alias	GALA; Galactosidase alpha (GLA)
Species	Human
Gene ID (Human)	2717
Gene ID	2717
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 150-429 of human Galactosidase alpha (Galactosidase alpha (GLA)) (NP_000160.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P06280
Protein Weight	49kDa
Shipping	Ice bag