

Anti-Galc Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P54803

Item Number: CM0027866

Introduction

High-quality anti-GALC rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Targets human, mouse, and rat galactocerebrosidase. Suitable for lysosomal storage disease research.

[Learn More](#)

Parameter	Value
Product Name	GALC Rabbit pAb
Remarks/Alias	GALC
Species	Human
Gene ID (human)	2581
Gene ID	2581
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-300 of human GALC (NP_000144.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P54803
Protein Weight	77kDa
Shipping	Ice bag