

Anti-Gale (Udp-Glucose 4-Epimerase) Polyclonal Antibody For Wb, Ihc-P, Elisa - Q14376

Item Number: CM0030391



Introduction

GALE Rabbit pAb validated for WB, IHC-P, ELISA. Detects human, mouse, rat GALE (UDP-glucose 4-epimerase). Ideal for galactose metabolism and glycoprotein research.

[Learn More](#)

Parameter	Value
Product Name	GALE Rabbit pAb
Remarks/Alias	SDR1E1; GALE
Species	Human
GeneID (Human)	2582
GeneID	2582
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 129-348 of human GALE (NP_001121093.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14376
Protein Weight	38kDa
Shipping	Ice bag