

Anti-Galectin-1 Polyclonal Antibody For Wb, If/Icc, Elisa - P09382

Item Number: CM0020165



Introduction

Rabbit polyclonal antibody targeting human Galectin-1 (LGALS1), validated for Western blot, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Useful for studying apoptosis, cell proliferation, and immune regulation.

[Learn More](#)

Parameter	Value
Product Name	Galectin 1/LGALS1 Rabbit pAb
Remarks/Alias	GBP; GAL1; Galectin 1/LGALS1
Species	Human
GeneID (Human)	3956
GeneID	3956
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-135 of human Galectin 1/LGALS1 (NP_002296.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09382
Protein Weight	15kDa
Shipping	Ice bag