

Anti-Galectin 1/Lgals1 Polyclonal Antibody For Wb, Elisa - P09382

Item Number: CM0020037



Introduction

This rabbit polyclonal antibody targets human galectin-1 (LGALS1), a key lectin involved in apoptosis and immune regulation. Validated for Western blotting (WB) and ELISA, it also cross-reacts with mouse and rat samples, supporting translational research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Galectin 1/LGALS1 Rabbit pAb
Remarks/Alias	GBP; GAL1; S1
Species	Human
GeneID	3956
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-135 of human Galectin 1/LGALS1 (NP_002296.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09382
Protein Weight	15kDa
Shipping	Ice bag