

Anti-Galectin 1/Lgals1 Rabbit Monoclonal Antibody For Wb, IHC-P, Elisa - P09382

Item Number: CM0002828



Introduction

High-quality rabbit monoclonal antibody against human Galectin 1/LGALS1, validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for apoptosis, immune regulation, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Galectin 1/LGALS1 Rabbit mAb
Remarks/Alias	GBP; GAL1; Galectin 1/LGALS1
Species	Human
GeneID (human)	3956
GeneID	3956
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-135 of human Galectin 1/LGALS1 (NP_002296.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09382
Protein Weight	15kDa
Shipping	Ice bag