


Anti-Galectin 3/Lgals3 Rabbit Monoclonal Antibody For Wb, If/Icc, Ip, Elisa - P17931

Item Number: CM0006739



Introduction

Rabbit monoclonal antibody targeting human galectin-3 (LGALS3) for WB, IF/ICC, IP, ELISA. Cross-reacts with mouse and rat. Recombinant full-length immunogen; detects 26 kDa protein. Ideal for inflammation, autophagy, and cancer studies.

[Learn More](#)

Parameter	Value
Product Name	Galectin 3/LGALS3 Rabbit mAb
Remarks/Alias	L31; GAL3; MAC2; CBP35; GALBP; GALIG; LGALS2; Galectin 3/LGALS3
Species	Human
GeneID (Human)	3958
GeneID	3958
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human Galectin 3/LGALS3 (NP_002297.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17931
Protein Weight	26kDa
Shipping	Ice bag