


Anti-Galectin 3/Lgals3 Rabbit Monoclonal Antibody For Wb, If/Icc, Ip, Fc, Elisa - P17931

Item Number: CM0008421



Introduction

Rabbit monoclonal antibody against human Galectin-3/LGALS3, validated for WB, IF/ICC, IP, FC, ELISA. Detects human and mouse Galectin-3. Suitable for studying lectin-mediated processes, inflammation, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Galectin 3/LGALS3 Rabbit mAb
Remarks/Alias	L31; GAL3; MAC2; CBP35; GALBP; GALIG; LGALS2
Species	Human
Gene ID (Human)	3958
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human Galectin 3/LGALS3 (NP_002297.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IP, FC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P17931
Protein Weight	26kDa
Shipping	Ice bag