


Anti-Galectin 4/Lgals4 Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P56470

Item Number: CM0007949



Introduction

Rabbit monoclonal antibody against human Galectin 4/LGALS4, validated for WB, IF-P, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for junctional assembly and adhesion research.

[Learn More](#)

Parameter	Value
Product Name	Galectin 4/LGALS4 Rabbit mAb
Remarks/Alias	GAL4; L36LBP; Galectin 4/LGALS4
Species	Human
Gene ID (Human)	3960
Gene ID	3960
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Galectin 4/LGALS4 (P56470)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P56470
Protein Weight	36kDa
Shipping	Ice bag