


Anti-Galectin 4/Lgals4 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P56470

Item Number: CM0021987



Introduction

Rabbit polyclonal antibody against human Galectin-4, validated for WB, IF-P, and ELISA. Cross-reacts with mouse and rat. Suitable for studying galactose-binding lectin functions and adherens junction biology.

[Learn More](#)

Parameter	Value
Product Name	Galectin 4/LGALS4 Rabbit pAb
Remarks/Alias	GAL4; L36LBP; Galectin 4/LGALS4
Species	Human
Gene ID (Human)	3960
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human Galectin 4/LGALS4 (NP_006140.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P56470
Protein Weight	36kDa
Shipping	Ice bag