

Anti-Galns Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P34059

Item Number: CM0026419



Introduction

High-quality rabbit polyclonal antibody against human N-acetylgalactosamine-6-sulfatase (GALNS), suitable for WB, IF/ICC, and ELISA. Detects human, mouse, and rat GALNS. Ideal for lysosomal storage disease research, including mucopolysaccharidosis IVA.

[Learn More](#)

Parameter	Value
Product Name	GALNS Rabbit pAb
Remarks/Alias	GAS; MPS4A; GaIN6S; GALNAC6S; GALNS
Species	Human
Gene ID	2588
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 353-522 of human GALNS (NP_000503.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P34059
Protein Weight	58kDa
Shipping	Ice bag