

Anti-Galnt2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q10471

Item Number: CM0027567



Introduction

High-quality rabbit polyclonal antibody targeting human GALNT2 (GalNAc-T2), validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for O-glycosylation research.

[Learn More](#)

Parameter	Value
Product Name	GALNT2 Rabbit pAb
Remarks/Alias	CDG2T; GalNAc-T2; GALNT2
Species	Human
Gene ID (Human)	2590
Gene ID	2590
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 442-571 of human GALNT2 (NP_004472.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q10471
Protein Weight	65kDa
Shipping	Ice bag