

Anti-Galnt3 Polyclonal Antibody For Wb And Elisa - Q14435

Item Number: CM0026055

Introduction

High-quality rabbit polyclonal antibody targeting human GALNT3 (GalNAc-T3), validated for Western blot (WB) and ELISA applications. Ideal for studies on O-linked glycosylation, mucin biology, and FGF23 regulation.

[Learn More](#)

Core Parameter	Specification
Product Name	GALNT3 Rabbit pAb
Remarks/Alias	HHS; HFTC; HFTC1; GalNAc-T3; GALNT3
Species	Human
Gene ID (Human)	2591
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human GALNT3 (NP_004473.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q14435
Protein Weight	73kDa
Shipping	Ice bag