

Anti-Galnt3 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q14435

Item Number: CM0018498

Introduction

Rabbit polyclonal antibody specific for human GALNT3 (GalNAc-T3). Validated in WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Ideal for O-glycosylation and FGF23-related disease research.

[Learn More](#)

Parameter	Value
Product Name	GALNT3 Rabbit pAb
Remarks/Alias	HHS; HFTC; HFTC1; GalNAc-T3; GALNT3
Species	Human
Gene ID (Human)	2591
Gene ID	2591
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 43-160 of human GALNT3 (NP_004473.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14435
Protein Weight	73kDa
Shipping	Ice bag