

Anti-Galnt6 Polyclonal Antibody For Wb, Elisa - Q8Ncl4

Item Number: CM0017579

Introduction

Rabbit polyclonal antibody against human GALNT6 (GalNAc-T6), validated for Western blot and ELISA. Designed for studying O-linked glycosylation, oncofetal fibronectin synthesis, and FGF23 modification.

[Learn More](#)

Parameter	Value
Product Name	GALNT6 Rabbit pAb
Remarks/Alias	GalNAcT6; GALNAC-T6; GALNT6
Species	Human
GeneID (Human)	11226
GeneID	11226
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 30-140 of human GALNT6 (NP_009141.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8NCL4
Protein Weight	71kDa
Shipping	Ice bag