

Anti-Galnt7 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q86Sf2

Item Number: CM0018997



Introduction

GALNT7 rabbit polyclonal antibody validated for WB, IHC-P, ELISA. Specifically detects human N-acetylgalactosaminyltransferase 7 (75kDa), a Golgi-resident glycopeptide transferase essential for O-linked oligosaccharide biosynthesis. Suitable for glycobiology and cancer studies.

[Learn More](#)

Parameter	Value
Product Name	GALNT7 Rabbit pAb
Remarks/Alias	GalNAcT7; GALNAC-T7; GALNT7
Species	Human
Gene ID (Human)	51809
Gene ID	51809
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-130 of human GALNT7 (NP_059119.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q86SF2
Protein Weight	75kDa
Shipping	Ice bag