

Anti-Galp Polyclonal Antibody For Wb, Elisa - Q9Ubc7

Item Number: CM0020628

Introduction

Rabbit polyclonal antibody against human GALP, suitable for Western blot and ELISA. Recognizes human, mouse, and rat GALP. Investigate neuropeptide function, hypothalamic signaling, and antimicrobial activity.

[Learn More](#)

Parameter	Value
Product Name	GALP Rabbit pAb
Remarks/Alias	GALP
Species	Human
Gene ID (human)	85569
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-116 of human GALP (NP_149097.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UBC7
Protein Weight	13kDa
Shipping	Ice bag