

Anti-Galt Polyclonal Antibody For Wb, Elisa - P07902

Item Number: CM0018920



Introduction

Rabbit polyclonal antibody targeting human GALT (galactose-1-phosphate uridylyltransferase). Suitable for WB and ELISA; cross-reacts with mouse and rat. Predicted molecular weight 43 kDa. Useful for galactose metabolism research.

[Learn More](#)

Parameter	Value
Product Name	GALT Rabbit pAb
Remarks/Alias	GALT; galactose-1-phosphate uridylyltransferase
Species	Human
Gene ID	2592
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-379 of human GALT (NP_000146.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07902
Protein Weight	43kDa
Shipping	Ice bag