

Anti-Gamma-Synuclein Monoclonal Antibody For Wb - O76070

Item Number: CM0011100

Introduction

Rabbit monoclonal antibody raised against human gamma-Synuclein (O76070). Validated for Western blot, IHC-P, and ELISA. Cross-reacts with mouse and rat. Suitable for neurofilament integrity and neurodegenerative disease research.

[Learn More](#)

Parameter	Value
Product Name	gamma-Synuclein Rabbit mAb
Remarks/Alias	SR, BCSG1, gamma-Synuclein
Species	Human
GeneID (Human)	6623
GeneID	6623
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 48-127 of human gamma-Synuclein (O76070).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O76070
Protein Weight	13kDa
Shipping	Ice bag