

Anti-Gamma-Synuclein Polyclonal Antibody For Wb, Ihc-P, Elisa - 076070

Item Number: CM0018509



Introduction

Rabbit polyclonal antibody targeting human gamma-Synuclein (SNCG/Breast cancer-specific gene 1). Validated for Western blot, IHC-P, and ELISA. Cross-reacts with human and mouse. Ideal for neurobiology and oncology research.

[Learn More](#)

Parameter	Value
Product Name	gamma-Synuclein Rabbit pAb
Remarks/Alias	SR; BCSG1; gamma-Synuclein
Species	Human
Gene ID	6623
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	076070
Protein Weight	13kDa
Shipping	Ice bag