

Anti-Ganab Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q14697

Item Number: CM0016262



Introduction

Rabbit polyclonal antibody against human GANAB, the catalytic subunit of glucosidase II. Validated for WB, IF/ICC, and ELISA; cross-reacts with mouse and rat. Ideal for glycoprotein processing and polycystin maturation studies.

[Learn More](#)

Parameter	Value
Product Name	GANAB Rabbit pAb
Remarks/Alias	G2AN; GIIA; PKD3; GLUII; GANAB
Species	Human
Gene ID (Human)	23193
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human GANAB (NP_001265123.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14697
Protein Weight	107kDa
Shipping	Ice bag