

Anti-Ganc Rabbit Polyclonal Antibody For Wb, Elisa - Q8Tet4

Item Number: CM0020824

Introduction

Rabbit polyclonal antibody against human neutral alpha-glucosidase C (GANC), validated for Western blot and ELISA. Shows cross-reactivity with mouse and rat. Suitable for glycosidase research.

[Learn More](#)

Parameter	Value
Product Name	GANC Rabbit pAb
Remarks/Alias	GANC
Species	Human
Gene ID (Human)	2595
Gene ID	2595
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human GANC (NP_937784.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q8TET4
Protein Weight	104kDa
Shipping	Ice bag