

Anti-Gap43 Monoclonal Antibody For Wb, If, Elisa - P17677

Item Number: CM0002459

Introduction

Rabbit monoclonal antibody against human GAP43 (Growth-associated protein 43, B-50, pp46) for Western blot, immunofluorescence, and ELISA; cross-reacts with mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	GAP43 Rabbit mAb
Remarks/Alias	B-50; PP46; GAP-43; GAP43
Species	Human
GeneID (Human)	2596
GeneID	2596
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 139-238 of human GAP43 (P17677).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-F, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17677
Protein Weight	25 kDa
Shipping	Ice bag