

Anti-Gap43 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P17677

Item Number: CM0028345

Introduction

Rabbit anti-human GAP43 polyclonal antibody validated for WB, IF/ICC, and ELISA. Detects GAP43/neuromodulin (25 kDa) across human, mouse, and rat. Ideal for neuroscience research on axon growth, synaptic plasticity, and neurodevelopment.

[Learn More](#)

Parameter	Value
Product Name	GAP43 Rabbit pAb
Remarks/Alias	B-50; PP46; GAP-43; GAP43
Species	Human
GeneID (human)	2596
GeneID	2596
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-238 of human GAP43 (NP_002036.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17677
Protein Weight	25kDa
Shipping	Ice bag