

Anti-Gapvd1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q14C86

Item Number: CM0021719

Introduction

GAPVD1 Rabbit pAb is a polyclonal antibody for WB, IF/ICC, ELISA with cross-reactivity to Mouse. GAPVD1 is a GTPase-activating protein and GEF involved in endocytosis and GLUT4 trafficking.

[Learn More](#)

Parameter	Value
Product Name	GAPVD1 Rabbit pAb
Remarks/Alias	RAP6; GAPEX5; GAPex-5; GAPVD1
Species	Human
Gene ID	26130
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human GAPVD1 (NP_056450.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse
SWISS	Q14C86
Protein Weight	165kDa
Shipping	Ice bag